Child Talker Beta How To Use:

Getting Started	Page 2
Types of Tiles	Page 3
Talker Tiles	Page 3
Folder Tiles	Page 3
Functional Tiles	Page 3

Getting Started

- I. Launch the application by clicking the "Child Talker Beta" Icon on the desktop
- II. The application will then launch and bring you to a screen with a red background
- III. To get started select "insert" in the upper left hand corner
 - A. If "Talker Tile" is selected user will be prompted to:
 - 1. Enter a phrase for the Talker Tile
 - 2. Select an image for the Talker Tile
 - 3. If both are provided a new Talker Tile will be added to the screen
 - B. If "Folder Tile" is selected user will be prompted to:
 - 1. Enter a name for the Folder Tile
 - 2. If a name is provided the folder tile will be added to the screen
- IV. To delete a tile right click on the Tile then select yes when asked if you would like to delete the Tile
- V. The application autosaves, any changes that are made will be saved on exit.

Types of Tiles

Talker Tiles (aka Items)

- I. Talker Tiles are Tiles that sound the phrases written on them when selected
- II. Talker Tiles can be deleted:
 - A. If a Talker Tile is right clicked then the user will be asked "Would you like to delete this tile?"
 - 1. If "Yes" is selected then the Tile will be deleted
 - 2. Otherwise if "No" is selected then the Tile will not be deleted.

Folder Tiles

- I. Folder Tiles are Tiles that contain another view / other Tiles
- II. When a Folder Tile is selected The Display will change to show the Tiles that were stored within the selected Folder.
 - A. When a Folder Tile is selected a "Back" Tile will be available to allow the user to return to the previous view.
- III. Folder Tiles can be deleted:
 - A. If a Folder Tile is right clicked then the user will be asked "Would you like to delete this tile?"
 - 1. If "Yes" is selected then the Tile will be deleted
 - 2. Otherwise if "No" is selected then the Tile will not be deleted.

Functional Tiles

- I. The "Back" Tile is a Functional Tile -> returns to the previous page
 - A. If you select "Back" you will return to the view you were at before selecting the most recent Folder Tile.
 - B. "Back" is only available when there is a view to go back to.